

**AIDMICS**  
BIOTECHNOLOGY

 iSperm

# User Manual

## iSperm 6

Version 6.3.0



# Index

<b>01 Hardware</b>	02
Specifications	03
Product Components	04 - 05
Installation (iSperm App, Optical Lens, Sample Collector, Heater)	06 - 14
Latex Beads Test	15 - 18
<b>02 Semen Dilution &amp; Sampling</b>	19
Basic Principle	20 - 23
Dilution Methods	24 - 25
Sampling Methods	26 - 33
<b>03 Software</b>	34
Semen Analysis	35 - 40
Semen Extending	41 - 45
Calculation Result	46 - 49
Supporting System	50 - 53
Data Center	54 - 56
Settings & Help	57 - 58
<b>04 Data Backup</b>	59
iSperm Cloud	60 - 61
iCloud	62 - 63
<b>05 Data Transfer</b>	64
Export Videos to Photo App	65 - 68
Export Data to Files App	69
Transfer Data from iPad	70 - 74
<b>06 Frequently Asked Questions</b>	75 - 81
<b>07 Validation</b>	82 - 88
<b>08 Certificate</b>	89 - 92

# 01

## Hardware



# Specifications


## Hardware

1. Optical magnification: equivalent to 200x in a traditional microscope.
2. Optical resolution: 1-1.5 $\mu$ m.
3. Heater: 37  $\pm$  0.5°C (DC 5V).
4. Weight: 350g (excluding the Apple iPad mini).
5. Power and battery: LR44 x 3, up to 45 hours.
6. Color: black.
7. Camera specs for iPad mini6.
  - 12MP photos.
  - 1080p/60 Fps Full HD recording.

## Software

Range of Analysis:

1. Concentration:



2. Motility :
  - 0%-100% for any concentration (**optimized at <500 million/ml**).
3. Progressive motility:
  - 0%-100% between concentration 10-75 million/ml (**optimized at 30-60 million/ml**).

Analysis Time: Concentration & Motility <20sec;

Progressive Motility ~30sec.

## Semen Sample

1. Fresh semen:
  - Direct measurement of raw semen for quick screening.
  - **Clear/Purified** extender (at 36-37°C) for diluted semen.
2. Thawed semen:
  - **Clear/Purified** extender (at 36-37°C) for frozen-thawed semen.

# Product Components



- 1 iSperm Briefcase
- 2 Sample Collector
- 3 Heater and Heater Cable
- 4 Hex Wrench
- 5 Sample Chips (Base & Cover)
- 6 Droppers
- 7 User Manual
- 8 Serial Number Card
- 9 iPad mini and iPad mini Case
- 10 Measuring Cups
- 11 LR44 Batteries
- 12 iPad Stand
- 13 Air Blower
- 14 Bag Strap

# Sampling Chips



15 500 Tests per box.

# Install iSperm App – 1/2

Search “iSperm 6” on App Store and download the species which needs to analysis.



iSperm Equine 6



iSperm Swine 6



iSperm Bovine 6



iSperm Caprine 6



iSperm Canine 6



iSperm Ovine 6




iSperm Poultry 6




iSperm Cervine 6

# Install iSperm App – 2/2


Open iSperm App. Follow the steps to **ACTIVATE** the serial number, provided in the Serial Number Card and **CREATE A NEW ACCOUNT or LOG IN** for iSperm Cloud.



# Install iSperm Case with Optical Lens




Insert the bottom corners into the case.




Then press the top corners (where LENS is located) of the iPad into the case.

**Reverse the order when removing the case.  
This prevents the case from being broken at the corner of the lens.**

# Use of Sample Collector




Press the silver ring and screw the end cap counterclockwise, and insert three new LR44 batteries with **the anode (+) facing up**. Screw back the end cap clockwise.



**If you see the light flashing instead of continuing on, please change the batteries.**


# Lens Alignment – 1/4

Adjustment 3-axis lens holders to align the lens for accurate measurements.






# Lens Alignment – 2/4



Turn on the light of the Sample Collector.

Screw the Sample Collector into the optical lens.




Open the iSperm App and tap “Lens Alignment.”


**Users are suggested to execute calibration before the first analysis.**

# Lens Alignment – 3/4

## Decentered lens



## Calibrated



Adjust the Optical Center + to the Target Center [ ] .

“Sample Collector is calibrated” shows up when the Optical Center overlaps the Target Center.


**Decentered Lens Detected**  
Decentered lens can lead to inaccurate analysis. Please adjust the lens to the center immediately.

Ignore Adjust Now


**Once entering the analyzing mode, the decentered lens will cause a warning message.**

# Lens Alignment – 4/4


**CAUTION: Use hex wrench gently**




Adjustment example 1:




**Original**




**Alignment approach**




Adjustment example 2:




**Original**




**Alignment approach**



Adjustment example 3:




**Original**



**Alignment approach**

# Heater Installation



1. Install the Heater with the optical lens holder.
2. Use a hex wrench to fasten the heater.
3. Connect the power cable to the heater.
4. Connect the power cable to a plug or a power bank.



**The LED indicator will turn on, once the heater is connected to the power. When the LED starts flashing, this means that the temperature has reached  $37\pm 0.5$  °C.**

# Latex Beads Test – 1/4

Latex Beads Test helps users to obtain more accurate readings with their iSperm. It will help to

- Check if the sampling process is correctly performed.
- Confirm if an abnormal hardware problem (sample collector, lens, ...) encounters.

For more information about latex beads, please refer to <https://www.hamiltonthorne.com/index.php/accu-beads>

## Preparation for Latex Beads Test

1. Familiarize yourself with the “Sampling > Pipette Method” in this manual to proceed with the test.
2. The **recommended concentration** interval for the latex beads solution to be prepared is 10 to 75 M/ml, and the size is 4  $\mu\text{m}$ .
3. 46 M/ml is the deal concentration for the test.  
(Refer to the picture on the right.)

Calibration by Chamber Type		
	For Fixed Coverslip	For Hemocytometer
#1	35 $\pm$ 5 M/ml	46 $\pm$ 7 M/ml
#2	18 $\pm$ 2.5 M/ml	23 $\pm$ 4 M/ml
#3	3 $\pm$ 1 M/ml	4 $\pm$ 1.5 M/ml

**Remind:**  
**Refer to the numbers shown on the column “For Hemocytometer” to do the comparison.**


## Instructions for Latex Beads Test

1. Make a sample and count the concentration with the “Latex Beads Test” Mode.
2. Count another aliquot of the beads sample with an iSperm Chip. The results should be within 10% of each other to be considered valid.
3. If the results are valid, average the two concentrations and compare them with the beads' acceptable ranges.
4. Record all results along with pertinent information, such as the person's name performing the procedure.


# Latex Beads Test – 2/4

## Instructions for Latex Beads Test

1. Open the iSperm app. Go to “Settings & Help” and tap “Latex Beads Test”.





2. Check that the view is OK, and tap “Start Analysis” to start the test.




# Latex Beads Test – 3/4

- Use 4-view analysis to obtain the averaged reading by rotating the arrow 90° for every measurement ( [See 4-view analysis detail in Software > Semen Analysis Section](#) ), then tap “Complete Analysis” to finish the tests.



# Latex Beads Test – 4/4

- Toggle to see the beads labeled in red, and the average of the four concentrations.



## Concentration too high or too low

Follow Step 4 to see if all the beads are well labeled.

### • If well-labeled

- Mixing error: Mix latex beads again for even distribution.
- Pipetting error: Reload the chip and be careful to avoid overloading the chip or underloading the chip.
- Confirm the latex beads' concentration using Hemocytometer.
- Contact the distributor if all the above issues are checked, and the problem still exists.

### • If not well-labeled







- It may be caused by lens contamination or abnormal light source. Contact the distributor for further assistance.




# 02

## Semen Dilution & Sampling

# Basic Principle - Preferred Specimens



# Basic Principle - Items in Use



**iSperm  
(iPad mini/Lens/Heater)**



Sample Collector



Stand




Base Chip




Eppendorf



Dropper



Pipette & Tips



Cover Chip











Measuring Cup

# Basic Principle - Semen Dilution

## Range of analysis:

Parameters	Full Range	Optimized Range
Concentration	10-1500 million/ml 500-1500 million/ml for quick preview test only	<500 million/ml
Motility	0%-100%	0%-100% at concentration <500 million/ml
Progressive Motility	0%-100% at concentration between 10 and 75 million/ml	0%-100% at concentration between 30 and 60 million/ml


Required diluting semen to 30-50 M/ml to obtain accurate readings and the kinematic parameters.

Species	Typical Concentration	Concentration that dilute (Semen : Extender, 1 : x)	
		Progressive test	Quickly test
 Equine	70-200 M/mL	1:5	1:1
 Canine	70-200 M/mL	1:5	1:1
 Swine	70-200 M/mL	1:5	1:1
 Bovine	800-1100 M/mL	1:30	1:5
 Caprine	1000-3000 M/mL	1:40	1:5
 Ovine	1000-3000 M/mL	1:40	1:5
 Cervine	1500-2500 M/mL	1:40	1:5
 Poultry	3000-5000 M/mL	Not Recommended	1:9

# Basic Principle - Dilution Ratio

## Dilution Ratio example:

Add semen into the extender is recommended.



Dilution Ratio	Semen	Extender
1:1 (2x)	100 µL	100 µL
1:2 (3x)		200 µL
1:3 (4x)		300 µL
1:4 (5x)		400 µL
1:9 (10x)	20 µL	180 µL
1:14 (15x)		280 µL
1:19 (20x)		380 µL
1:29 (30x)		580 µL
1:39 (40x)		780 µL

# Dilution Methods – 1/2


## • Partial Dilution (For Progressive Motility)

Taking a small portion of semen and diluting it with an extender. Dilution ratio ranges from 1:5 to 1:10.

Commonly used when :

1. Progressive motility readings are needed.
2. Raw semen is limited and needs to be preserved for AI (e.g., Equine, Poultry, Bovine, Ovine).

Tools Needed: Eppendorf, Micropipette (10-1000 $\mu$ L).



# Dilution Methods – 2/2

- **All-in Dilution (Convenient, Fast)**


Dilution of the entire collected raw semen with an extender.

Dilution ratio used is 1:1 or 1:2.

Commonly used when :

1. An extender is added to prolong sperm motility.
2. Raw semen is sufficient, meeting AI standards even after dilution (e.g., Swine, Canine).


Tools Needed: Beaker.



# Preparatory Work


**This step is necessary for all sampling methods.**

- Use the working board in the iSperm Briefcase, or find a flat, firm/stable surface (e.g., table).
- Clean the working board/table surface.
- Dusty surface could contaminate “Cover Chip” and hinder the analysis.



**iSperm view of dust/  
fiber on the Cover Chip.**

- Mount “Base Chip” onto Sample Collector.
- Place “Sample Collector” on the table.






# Three Methods of Sampling


## 1 Pipette Method

- Sample volume: 7.5 $\mu$ L.
- Best when the technician is familiar with pipette skills or when limited sperm is available.




## 2 Dropper Method

- Sample volume:  $\approx$ 70 $\mu$ L.
- Best when preferring not to handle pipette and with sufficient sperm.




## 3 Dipping Method

- For quick preview tests or for fixed sperm.




# Pipette Method – 1/2

1. Use a Pipette to mix the semen gently and thoroughly.
  - **10-15 times** to make the semen well-mixed and evenly distributed.
  - **Gentle speed** to prevent bubbles.
2. Drop **7.5µL** onto the Base Chip center area.



3. Give Sample Collector a gentle knock against the table to spread the droplet.

- **Flat semen droplet is crucial to reduce CV.**




a gentle knock would do the trick

**\*Go to next step as soon as possible**


## Pipette Method – 2/2

4. Place **“Cover Chip”** on a clean table with open-side facing up.
5. Flip Sample Collector downward.
  - **Droplet will remain on Base Chip.**
6. Press Base Chip into Cover Chip vertically. One will hear a **“click”** first; then, continue to press down for another 1-2 seconds.




# Dropper Method – 1/2

1. Use a dropper to mix the semen gently and thoroughly.
  - **10-15 times** in general to make the semen well-mixed and evenly distributed.
  - **gentle speed** to prevent bubbles.




2. Use a dropper to inject **70-100 $\mu$ L** into Cover Chip.
  - About **2mm height** (~100 $\mu$ L) in Cover Chip is sufficient.
  - **Over 3mm height** (~150 $\mu$ L) will lead to spillover when enclosing with Base Chip.



# Dropper Method – 2/2


3. Place semen-filled “Cover Chip” on a clean table with open-side facing up.
4. Sample Collector Flip downward.
5. Press Base Chip into Cover Chip vertically. One will hear a “click” first; then, continue to press down for another 1-2 seconds.



# Dipping Method – 1/2


This method is for quick preview tests or fixed semen.

1. Use a high volume Dropper/Pipette to mix the semen gently and thoroughly.
2. Wait 30 seconds to let the flow settle.
  - **Fast** microflow causes uncertainty during the dip sampling.




3. Dip Base Chip into semen.

- Immerse Base Chip **5mm below the semen surface**.



## Dipping Method – 2/2

4. Place **“Cover Chip”** on a clean table with open-side facing up.
5. Flip Sample Collector downward; **Droplet will remain on Base Chip.**
6. Press Base Chip into Cover Chip vertically. One will hear a **“click”** first; then, continue to press down for another 1-2 seconds.



# 03

## Software



# Semen Analysis – 1/4


1. Place the iPad on a stand.



**Do NOT lean Sample Collector against the table**


◀ Stand

2. Open the iSperm app, and tap “Analyze Now.”




# Semen Analysis – 2/4

3. Preview image. ( see details in “Preferred Specimens” )





# Semen Analysis – 3/4

4. Rotate Collector until the logo is pointing up. Tap “Start Analysis.”  
**Multiple-view (typically 4-view) analysis is recommended.**



**Rotating  
Collector gently  
will do the trick**

5. Check the analyzed result.



**Analyzed  
Results**


**Tap “Show All Data”  
to view the details.**

**Tap “X” to delete  
a measurement.**


**Tap image to  
view a video**

# Semen Analysis – 4/4

- Rotate Collector to complete 4-view analysis. Tap “Start Analysis” for every direction.



- Check the “average” of the repeated tests on the screen. Tap “Complete Analysis” to finish and conclude the tests.



**The relationship between views and coefficient of variation (CV) values:**


	Quick preview test			Recommend
view	1	2	3	4
CV	<div style="display: flex; align-items: center; justify-content: center;"> <span style="margin-right: 10px;">higher CV</span> <span style="font-size: 2em;">←</span> <span style="margin-left: 10px;">▶</span> <span style="margin-left: 10px;">lower CV</span> </div>			

# Manual Assessment of Abnormal Morphology

During the “Semen Analysis” step, tap the icon to zoom in to check sperm morphology manually.

## Follow the steps:

1. Use the manual focus bar to adjust the focus.
2. Tap the “Pause button” to freeze the image.
3. Fill in the estimated abnormality into the text field.
4. Tap “Save” to leave.



# Average of Analysis Date

Tap once to stop the video.

Tap twice to view the fullscreen video.

Slide left to view other indices.

Tap to Draw / Hide Sperm Tracks

Sperm tracks are labeled in 4 colors:

● Progressively Motile

● Motile ● Static ● Late Track

The screenshot shows the AIMICS 01 app interface. At the top, there is a 'Cancel' button (X) and a 'Save' button (checkmark). The main area displays a video of sperm tracks with a 'Multiple Issues' warning. A 'Draw Sperm Tracks' toggle is visible. Below the video is a summary card titled 'Average of Analysis Data (4)' with the following data:

Concentration	Motility	Progressive	Abnormality
52.08 M/ml	80 %	74 %	0 %

Below this is a 'Sperm Kinetics' section with the following data:

VCL	VAP	VSL	STR	LIN	BCF	ALH	WOB
132 $\mu\text{m/s}$	95 $\mu\text{m/s}$	78 $\mu\text{m/s}$	81 %	62 %	14 Hz	13 $\mu\text{m}$	74 %

At the bottom, there is a 'Was the Semen Diluted?' section with three options: 'Partial Dilution', 'All-in Dilution' (selected), and 'Raw / Chilled / Frozen'. The 'All-in Dilution' option shows a 'Raw Concentration' of 104.16 million/ml. A 'Dilution Ratio, Semen : Extender = 1 : 1' is also displayed.

Annotations on the screenshot include: 'Leave without saving.' pointing to the Cancel button, 'Save the result.' pointing to the Save button, and a hand icon pointing to a blue sperm track in the video.

❗ **Severe impurity. Unable to track sperms correctly.**

⚠️ **Large bubbles or debris detected.**

⚠️ **Shake detected.**

Concentration, motility, and progressive motility are on the first row.

Detailed kinetics are on the second row.

# Semen Extending – 1/5

## Step. 1 Was the Semen Diluted?

### 1 Partial Dilution (Fill in dilution ratio)

Taking a small portion of semen and diluting it with an extender for testing. In most cases, the dilution ratio ranges from 1:5 to 1:10.

Raw Concentration = Measured Concentration \* Dilution Ratio

### 2 All-in Dilution (Fill in dilution ratio)

Dilution of the entire collected raw semen with an extender. In most cases, the dilution ratio is 1:1 or 1:2.

Raw Concentration = Measured Concentration \* Dilution Ratio

In this case, the dilution ratio is 1:1, with a Measured Concentration of 52.08 M/ml.

Therefore, the Raw Concentration will be 104.16 M/ml.

e.g.:  $104.16 = 52.08 \times (1+1)$

### 3 Raw/Chilled/Frozen

No further dilution is conducted before analysis.

Concentration = Measured Concentration

The screenshot shows a mobile application interface for semen analysis. At the top, there are four main metrics: Concentration (52.08 M/ml), Motility (80%), Progressive (74%), and Abnormality (0%). Below these is a 'Sperm Kinetics' section with various parameters: VCL (132 μm/s), VAP (95 μm/s), VSL (78 μm/s), STR (81%), LIN (62%), BCF (14 Hz), ALH (13 μm), and WOB (74%). The main section is titled 'Was the Semen Diluted?' and contains three options: 1. Partial Dilution (unchecked), 2. All-in Dilution (checked), and 3. Raw / Chilled / Frozen (unchecked). Each option has a 'Raw Concentration' field. For 'All-in Dilution', the value is 104.16 million/ml. Below the options is a 'Dilution Ratio, Semen : Extender = 1 : 1' field. At the bottom, there are two tabs: 'Analysis Result' and 'Semen Extending'. A blue box highlights the 'Semen Extending' tab, and a blue line points from the text below to it.

Concentration	Motility	Progressive	Abnormality
52.08 M/ml	80 %	74 %	0 %

VCL	VAP	VSL	STR	LIN	BCF	ALH	WOB
132 μm/s	95 μm/s	78 μm/s	81 %	62 %	14 Hz	13 μm	74 %

Option	Raw Concentration
1 Partial Dilution	million/ml
2 All-in Dilution	104.16 million/ml
3 Raw / Chilled / Frozen	million/ml

Dilution Ratio, Semen : Extender = 1 : 1

Analysis Result | Semen Extending

Tap “Semen Extending” or slide left to go to the Extending page.

# Semen Extending – 2/5

The choice made in **Step. 1** will influence the concentration display of “Analysis Result”.

If you choose:

## 1 Partial Dilution

Concentration =  
Raw Concentration or  
pre-analysis diluted  
concentration

In the case, for instance,  
the dilution ratio is 1:5, the  
concentration of “Analysis  
Result” will be 312.48 M/ml.  
e.g.:  $312.48 = 52.08 \times (1+5)$

## 2 All-in Dilution

Concentration =  
Measured Concentration

## 3 Raw/Chilled/Frozen

Concentration =  
Measured Concentration

Motility, Progressive  
and Abnormality  
remains the same.

The image shows two screenshots of the 'Semen Extending' software interface. The top screenshot is the 'Was the Semen Diluted?' selection screen, which has three options: 1. Partial Dilution (unchecked), 2. All-in Dilution (checked), and 3. Raw / Chilled / Frozen (unchecked). The 'All-in Dilution' option shows a raw concentration of 104.16 million/ml and a dilution ratio of 1:1. A blue arrow points from this screen to the bottom screenshot. The bottom screenshot is the 'Analysis Result' screen for 'AIDMICS 01'. It displays a concentration of 52.08 M/ml, motility of 80%, progressive motility of 74%, and abnormality of 0%. Below the results, there are sections for '1. Semen Volume for Dilution' (50 ml), 'Semen Summary' (Total Sperm = 2.6 billion, Motile Sperm = 2.08 billion, Progressively Motile Sperm = 1.94 billion), '2. Volume per Dose' (20 ml, Variable checked), and '3. Extend based on:' (Effective Sperm per Dose selected, 400 million / dose). A blue arrow points from the 'Motility, Progressive and Abnormality' text on the left to the corresponding values in the 'Analysis Result' screen.



# Semen Extending – 3/5

## Step. 2-1 Semen Volume for Dilution

If you choose:

### 1 Partial Dilution

Please fill in the remaining volume available for extending after deducting the portion used for analysis.

### 2 All-in Dilution

Please fill the total volume of the semen and the added extender.

### 3 Raw/Chilled/Frozen

Please fill in the volume received.

## Semen Summary

This section displays how many sperm are available for extending :

### • Total Sperm

Total Sperm = Concentration x Semen Volume for Dilution  
e.g.: 2.6 billion = 52.08 million x 50 ml

### • Motile Sperm

Motile Sperm = Total Sperm x Motility(%)  
e.g.: 2.08 billion = 2.6 billion x 80%

### • Progressively Motile Sperm

Progressively Motile Sperm = Total Sperm x Progressive (%)  
e.g.: 1.94 billion = 2.6 billion x 74%

**Analysis Result**

Concentration	Motility	Progressive	Abnormality
52.08 M/ml	80 %	74 %	0 %

**1. Semen Volume for Dilution**

50 ml

**Semen Summary**

- Total Sperm = 2.6 billion
- Motile Sperm = 2.08 billion
- Progressively Motile Sperm = 1.94 billion

**2. Volume per Dose**

20 ml  Variable  Fixed

**3. Extend based on:**

Effective Sperm per Dose (million / dose)  Concentration after Dilution (M / ml)

400 million / dose  Total Sperm  Motile Sperm  Progressively Motile Sperm

Analysis Result | Semen Extending

# Semen Extending – 4/5

## Step. 2-2 Volume per Dose

- **Variable**

The container size is flexible (e.g., beakers, syringes).  
Allowing iSperm to adjust the dose volume to meet extending standards ( set on the **Step. 3** ) when semen quality is poor.

- **Fixed**

The container size is fixed (e.g. straws)  
iSperm is not allowed to change the dose volume.

**When semen quality is good, there is no difference between the two options.**

The screenshot shows the 'New Record' screen in the iSperm app. At the top, there are navigation options: 'Cancel', 'Support', and 'Save'. The date and time are '2023/12/13 18:17' and the device is 'EQ 6.3.12 (ADT7) | mini 6'. The main section is titled 'Analysis Result' and displays four key metrics: Concentration (52.08 M/ml), Motility (80 %), Progressive (74 %), and Abnormality (0 %). Below this, there are three sections: 1. Semen Volume for Dilution (50 ml), 2. Volume per Dose (20 ml), and 3. Extend based on. In section 2, the 'Variable' option is selected. In section 3, 'Effective Sperm per Dose' is selected, and the 'M / ml' unit is chosen. The 'Extend based on' section also includes checkboxes for 'Concentration after Dilution', 'Total Sperm', 'Motile Sperm' (which is checked), and 'Progressively Motile Sperm'. At the bottom, there are two tabs: 'Analysis Result' and 'Semen Extending'.

# Semen Extending – 5/5

## Step. 3 Extend based on

iSperm must follow this standard when extending to ensure the quality of each dose.

### • Effective Sperm per Dose

Type-in the required quantity and select the sperm type

In this case, 400 million sperm with motility per dose is needed.

\*When there's no progressive motility readings, the option will not be available for selection.

**3. Extend based on:**

Effective Sperm per Dose  Concentration after Dilution

million / dose M / ml

400 million / dose

Total Sperm  
 Motile Sperm  
 Progressively Motile Sperm

**Calculation Result**

Volume for Extender	Total Extender Volume	Number of Doses
50 ml	+ 61.6 ml	5
Semen Per Dose	Extender per Dose	
7.68 ml	+ 12.32 ml	

Summary per Dose

Choose Total, Motile or Progressive Sperms for calculation.

### • Concentration after Dilution

Enter the desired concentration for each dose

This value should be lower than the concentration in the “Analysis Result” on top, or the semen sample won’t be extendable. (in cases of low semen concentration)

**3. Extend based on:**

Effective Sperm per Dose  Concentration after Dilution

million / dose M / ml

0.75 M/ml

**Calculation Result**

Volume for Extender	Total Extender Volume	Number of Doses
50 ml	+ 2.07 ml	69
Semen Per Dose	Extender per Dose	
0.72 ml	+ 0.03 ml	

Summary per Dose

# Calculation Result – 1/4

## Calculation Result

- **Volume for Extension**

Volume entered in “1. Semen Volume for Dilution”

- **Total Extender Volume**

The total amount of extender required for extension

Total Extender Volume = Extender Per Dose x Number of Doses

e.g.:  $61.6 = 12.32 \times 5$

- **Number of Doses**

Number of Doses = Total or Motile or Progressively Motile Sperm / Effective Sperm per Dose

e.g.:  $5 = 2.08 \text{ billion} / 400 \text{ million}$

- **Semen Per Dose**

Semen Per Dose = Effective Sperm per Dose / Concentration

e.g.:  $7.68 = 400 / 52.08$

- **Extender Per Dose**

Extender Per Dose = Volume per Dose - Semen Per Dose

e.g.:  $12.32 = 20 - 7.68$

The screenshot shows a mobile application interface for a "New Record" entry. At the top, it displays the date and time (1:13 PM, Thu Dec 14) and battery status (29%). The main form includes several input fields and checkboxes. A blue line points from the "Calculation Result" section to the "Effective Sperm per Dose" field.

Field	Value	Unit
Effective Sperm per Dose	400	million / dose
Concentration after Dilution	M / ml	
Total Sperm	<input type="checkbox"/>	
Motile Sperm	<input checked="" type="checkbox"/>	
Progressively Motile Sperm	<input type="checkbox"/>	

Calculation Result		
Volume for Extension	+	Total Extender Volume
50 ml		61.6 ml
Semen Per Dose	+	Extender per Dose
7.68 ml		12.32 ml
		Number of Doses
		5

Summary per Dose	
Volume	= 20 ml
Concentration	= 20 M / ml

# Calculation Result – 2/4

The section's color will indicate the quality standard:

**Green:** Meets the standard (as per the value entered in “3. Extend based on”)

**Red:** Does not meet the standard

## Case 1: Number of Doses > 1 (Able to be divided into multiple doses)

If selected :

- **Variable:**

Can be split into doses, be **Green**.

- **Fixed:**

If Semen Per Dose ≤ Volume per Dose, it will be **Green**; otherwise, it will be **Red**.

20 ml  Variable  Fixed

3. Extend based on:

Effective Sperm per Dose: 1200 million / dose

Concentration after Dilution: M/ml

Calculation Result:

Volume for Extension: 50 ml + Total Extender Volume: 0 ml = Percentage of Doses: 1

Semen Per Dose: 23.04 ml Extender per Dose: 0 ml

20 ml  Variable  Fixed

3. Extend based on:

Effective Sperm per Dose: 1200 million / dose

Concentration after Dilution: M/ml

Calculation Result:

Volume for Extension: 50 ml + Total Extender Volume: 0 ml = Number of Doses: 1

Semen Per Dose: 20 ml Extender per Dose: 0 ml

## Case 2: Number of Doses < 1 (Unable to obtain even one dose)

If selected:

- **Variable**

(All packed as one dose)

Number of Doses = 1

- **Fixed**

(Divide evenly)

$$\text{Number of Doses} = \frac{\text{Volume for Extension}}{\text{Volume per Dose}}$$

50 ml + 0 ml = 1


50 ml + 0 ml

50 ml + 0 ml = 2

20 ml + 0 ml

# Calculation Result – 3/4

## Flowchart




# Calculation Result – 4/4

## Summary Per Dose

### ● Volume: Volume of each dose

In Section “3. Extend based on”, if “Fixed” is selected, the dose volume should match the entered volume.

If “Variable” is selected and iSperm adjusts the volume, an **⚠ exclamation mark** appears. In such cases, please verify whether the readings for each parameter align with Section “3. Extend based on” standards.



**Summary per Dose**

Volume = **23.04 ml** ⚠

Concentration = **52.08 M / ml**

Total Sperm = **1041.63 million**

Motile Sperm = **1041.63 million**

Progressively Motile Sperm = **775.52 million**


The screenshot shows a 'Summary per Dose' card with a list of metrics. A blue callout box highlights the 'Volume' value of 23.04 ml, which is accompanied by a yellow exclamation mark icon.

### ● Concentration

### ● Total Sperm

### ● Motile Sperm

### ● Progressively Motile Sperm




The screenshot displays two sections of a mobile application interface. The top section, titled 'Summary per Dose', shows a list of metrics: Volume = 20 ml, Concentration = 20 M / ml, Total Sperm = 400 million, Motile Sperm = 320 million, and Progressively Motile Sperm = 297.81 million. A blue callout line points from the 'Volume' value in this section to the 'Volume' value in the callout box above. The bottom section, titled 'Information', contains input fields for 'Name' and 'Source', and a 'Notes' field with a placeholder text: '(This is a read-only field for displaying the notes of this analysis result)'. A blue callout line points from the 'Notes' field to the text below.

Fill in Name, Source and the supplemental information.

# Supporting System – 1/4

Tap Support to send the data of abnormal measurements via email directly from the iPad. (Please set up the Mail app before using this function.)



## Set Up the Mail Apps

If you're been using Mail App on your iPad, please skip to "Supporting System (3/4)".

1. Go to "Settings" and Scroll down the sidebar to find the "Mail" and tap "Accounts".






## Supporting System – 2/4

2. Tap “Add Account” and choose the service that you own an account on, such as Google or Yahoo.


**If this is a shared device, we suggest creating a new account for this device.**

**You may follow the guide below to create an iCloud Account:**

<https://support.apple.com/guide/icloud/create-an-icloudcom-email-address-mmdd8d1c5c/1.0/icloud/1.0>




4. After login, please turn on the “Mail” service with this account and click to save.



# Supporting System – 3/4

1. Tap “Support,” and the required data will be appended to an email draft automatically.



**If you see the pop-up on the right, go back to check if the mail sender is ready.**

## Sender Not Configured


The sender is not configured on this iPad. If you aren't sure how to set up an email address, please refer to iSperm software manual, or contact your local distributor.

OK



# Supporting System – 4/4

- Please describe the issue and **DO NOT DELETE** any auto-generated information.
- Make sure the Internet connection is okay, and tap to send the supporting mails.

- The “Email Sent” will pop up after the mail has been sent.




- After a few minutes, you will receive an auto generated support request mail in the Mail App. Please wait for further technical support.



**If you didn't see the mail, please check your Internet access and report again, or contact your distributor. (Auto-generated emails might be categorized as spam, so please check the spam mailbox as well.)**

# Data Center - Local Experiments



A list of measurement results stored on the local iPad.

Tap any measurement to view detailed results.


- Videos can be viewed if they are stored during the analysis.
- Only the Animal Name and Source can be edited.

Analysis ID	Concentration	Motility	Progressive	Status
AIDMICS 01	115.45 M / ml	85%	N/A	Unsynced
AIDMICS 02	0.22 M / ml	100%	0%	Unsynced
AIDMICS 03	0.00 M / ml	0%	0%	Video Synced
AIDMICS 04	3.73 M / ml	100%	81%	Synced
AIDMICS 05	161.48 M / ml	100%	N/A	Synced
AIDMICS 06	153.96 M / ml	100%	N/A	Synced
AIDMICS 07	103.67 M / ml	80%	N/A	Synced

Swipe to the left and tap “Delete” to delete unwanted measurement result.




# Data Center - Online Charts




A graph helps you to visualize all the measured history.

Tap to change X and Y-Axis.



Tap to change Sort.



Animal ID	Measure Time	Raw Conc.	Mottle	Progress
AIDMICS 01	2022-05-05 10:00:18	371.16 M/ml	92.79 M/ml	71.42 %
AIDMICS 02	2022-04-14 07:41:18	0 M/ml	0 M/ml	0 %
AIDMICS 03	2022-04-14 16:49:38			


Tap to download Measurement history in Excel format.

Tap to view the detailed results.

Abnormality	Dilution Ratio	Total Sperm	Effective Sperm	Total Volume
0 %	1:3	22.27 Billion	19.91 Billion	85 ml


Tap Animal Report to view all measurements.

# Data Center - Animal List



View measurement history of a specific ID.

Tap to change Sort.





Tap Animal Report to view all measurements.



Data shows the average from all measurements.

# Settings & Help – 1/2



The features of Settings are on the next page.

Email & Phone.

Log in/out iSperm Cloud account.

Owner info; Serial Number;

Information about system and app.

# Settings & Help – 2/2

## Analysis Settings:

- Round Number of Doses
- Motility Parameters

1. Adjust cutoff values of progressive sperm assessment.
2. Adjust cutoff values of motile sperm assessment.
3. Calculation Time

## Definitions



VCL	Curvilinear Velocity (µm/s)
VAP	Average Path Velocity (µm/s)
VSL	Straight Line Velocity (µm/s)
STR	Straightness (VSL/VAP)
LIN	Linearity (VSL/VCL)

**Round Number of Doses** ☑  
Should (Sperm) round number of doses to the nearest integer?

---

**Motility Parameters** ✎  
Progressive Cutoff Value: VAP = 116, STR = 66%  
Motility Cutoff Value: VAP = 50, VSL = 13  
Calculation Time: 45 frames

---

**Motility Parameters**

**1** **Progressive Cutoff Value**  
(Default: VAP = 60, STR = 75%)  
VAP = 116 µm/s (0-150)  
STR = 66 % (0-100)  
A sperm matching these conditions will be labeled progressive.

**2** **Motile Cutoff Value**  
(Default: VAP = 20, VSL = 0)  
VAP = 50 µm/s (0-116)  
VSL = 13 µm/s (0-50)  
Sperms matching these conditions will be labeled motile, and will be included in the calculation of the average VCL, VAP, VSL, STR and LIN.

**3** **Calculation Time**  
(Default: 35 frames, at 50 fps)  
Calculate with: 45 frames (10-90)  
This setting will affect motility, progressive motility and all velocity.

**VCL** (Curvilinear velocity) =  $\frac{\text{The summation of distance between the sperm head positions in each frame}}{\text{Elapsed time}}$

**VAP** (Average path velocity) =  $\frac{\text{The averaged path determined by smoothing the position of head}}{\text{Elapsed time}}$

**VSL** (Straight line velocity) =  $\frac{\text{Distance between the first and last points}}{\text{Elapsed time}}$

**STR** (Straightness) =  $\frac{\text{VSL/VAP in percent(\%)}}{\text{is a measure of track compactness.}}$

**LIN** (Linearity) =  $\frac{\text{VSL/VCL in percent(\%)}}{\text{is a measure of track direction.}}$



# 04

## Data Backup

- **iSperm Cloud**

iSperm Cloud allows users to view measurement history using web browsers (Chrome, Safari, etc.) on desktop computers, laptops, and cellular phones.


- **iCloud**

iSperm App supports a full backup for iPad on iCloud.

All the videos and data in the iSperm App can be restored if you need to use a new iPad.


# iSperm Cloud – 1/2

1. Login iSperm account at <https://ispermcloud.aidmics.com/>



The image shows the login page for iSperm Cloud. At the top center is the iSperm Cloud logo, which consists of a stylized 'S' icon followed by the text 'iSperm Cloud'. Below the logo are two input fields: 'Account (email)' with the placeholder text 'account@email.com' and 'Password' with masked characters '\*\*\*\*\*'. To the right of the password field is a 'Show' button with an eye icon. Below the input fields is a 'Login' button, followed by links for 'New Account' and 'Forgot Password?'.


2. Choose animal species.



The image shows the species selection page in iSperm Cloud. At the top left is the text 'iSperm Cloud' and at the top right are links for 'Animal List', 'Experiment List', and 'Your Account'. The main content area has the heading 'What species do you want to see?' followed by a row of ten circular icons representing different animal species: Bovine, Caprine, Cervine, Canine, Equine, Murine, Ovine, Poultry, and Swine. Each icon is labeled with its corresponding species name. At the bottom of the page, there is a navigation bar with a language selector showing 'en' (English) and '正' (Chinese), and a selected species icon for 'Equine'.

# iSperm Cloud – 2/2

3. All the functions and user interface on the “Data Center” are almost identical to the iSperm App.




# iCloud Backup – 1/2

**Before backup, be sure that you have enough space available in iCloud.**

When you sign up for iCloud, you will get 5GB of iCloud storage for free. If you need more iCloud storage, you can buy more from your iPad with your Apple ID. Learn more about prices in your region: <https://support.apple.com/kb/ht201238>

1. Connect your device to a Wi-Fi network.
2. Go to Settings > [your name], and tap iCloud.
3. Tap iCloud Backup.
4. Tap Back Up Now. Stay connected to your Wi-Fi network until the process completes.



# iCloud Backup – 2/2


## Check the status of the iCloud backup

You can check the progress and confirm if the backup is completed. Go to Settings > [your name] > iCloud > iCloud Backup. Under Back Up Now, you'll see the date and time of your last backup.

## Automatically back up with iCloud Backup

To let iCloud automatically back up your device each day, here's what you need to do:

- Make sure that iCloud Backup is turned on in Settings > [your name] > iCloud > iCloud Backup.
- Connect your device to a power source.
- Connect your device to a Wi-Fi network.
- Make sure that your device screen is locked.



# 05

## Data Transfer

- **Export Videos to Photo App**


Videos taken on the iSperm app can be exported to Photos App. Each video would come with its respective reading set (concentration, motility, ...) shown underneath the video.

- **Transfer Data from iPad**

You can copy videos and measurements between your computer and apps on your iOS device using File Sharing.

# Export Videos to Photo App – 1/4

1. Go to Local Experiments.  
Tap the Export Button to enter the “Selecting” Mode.



Time	Concentration	Mobility	Progressive	Status
2022/05/10 11:31	115.45 M/ml	85%	N/A	Unsynced
2022/05/10 11:39	0.22 M/ml	100%	0%	Unsynced
2022/05/10 11:59	0.00 M/ml	0%	0%	Value Synced
2022/05/10 11:56	3.73 M/ml	100%	81%	Synced
2022/05/10 10:58	161.48 M/ml	100%	N/A	Synced
2022/05/10 10:59	153.96 M/ml	100%	N/A	Synced
2022/05/10 10:29	103.67 M/ml	80%	N/A	Synced

Time	Concentration	Mobility	Progressive
2022/05/10 11:31	115.45 M/ml	85%	N/A
2022/05/10 11:39	0.22 M/ml	100%	0%
2022/05/10 11:58	0.00 M/ml	0%	0%
2022/05/10 11:56	3.73 M/ml	100%	81%
2022/05/10 10:58	161.48 M/ml	100%	N/A
2022/05/10 10:59	153.96 M/ml	100%	N/A
2022/05/10 10:29	103.67 M/ml	80%	N/A

## Data Status


**Unsynced** Unsynced data.

**Value Synced** Readings are synced, syncing the videos.

**Synced** Reading and videos are synced.

# Export Videos to Photo App – 2/4



2. Select the measurements you want to output.
3. Check the box if tracks are needed.  
( To see the difference, please refer to “Export Videos to Photo App - 4/4” )
4. Tap the Export button.






# Export Videos to Photo App – 3/4

5. Await until all videos are exported.




Export result



6. Go to Photos App.


Tap the iSperm's Albums and find the videos (with uneditable readings) to play.




# Export Videos to Photo App – 4/4

## Videos exported with/without Sperm Tracks

Export Sperm Tracks




Export Sperm Tracks




# Export Data to Files App

1. Check the box to export the video. (Follow the “Export videos to Photo App” step.2~5.)




2. Go to Files App. Tap On My iPad and find the videos and the editable readings.




# Transfer Data from iPad (iTunes) – 1/3

Please refer to the “Transfer Data From iPad (Finder)” if your computer is MacOS Catalina or later.

1. Open iTunes on your Mac or PC.
2. Connect your iPad to your computer using the USB cable that comes with your device. You might see a prompt on the iOS device asking you to Trust This Computer. Tap Trust to continue.




3. Click your device on iTunes.




# Transfer Data from iPad (iTunes) – 2/3

4. In the left sidebar, click File Sharing.
5. Select the iSperm App. Videos are stored in folders for each measurement. The folder name format is **“Name + Species + Date”**. You can save the videos to your Mac or PC.



# Transfer Data from iPad (iTunes) – 3/3

- When opening the folder from your Mac or PC, you will see multiple videos (depending on how many views you analyzed) and a “tsv” file that summarizes all the measurements.



NAME1\_2019-06-12\_08-52-17.tsv

Equipment	Arm Name	Programs Count	Matrix Count	Device Type	Time
NAME1		800 = 70 & VSP = 80	VSP = 20 & VSL = 0	N/A	2019-06-12 08:52:17

Address	ID	Concentration (Infectious)	Matrix (%)	Progressive (%)	VSP (uV/s)	VSL (uV/s)	VSL (uV/s)	LRN (%)	SIR (%)	Abnormal Tracks Detected	Shells Detected	Subsides / Large Debris Detected	Impurity
1	239.63	95	87	118.43	238.32	35.42	33.9	44.59	NO	YES	NO	Detected	
2	148.6	95	33	78.53	148.86	43.52	80.89	44.42	NO	YES	NO	Detected	
3	187.87	94	11	68.88	138.32	34.78	85.42	40.39	NO	YES	NO	Detected	
4	88.76	95	4	83.71	182.72	71.98	84.24	40.39	NO	YES	NO	Detected	
Median	239.29	94	22	82.22	188.79	47.88	82.2	43.3					

Packaging	Number of Droplets	Estimate Volume (nL)	Number of Spores (Billions)	Total Volume (nL)	Number of Matrix Spores per Droplet (Infectious)	Volume per Droplet (nL)	Dilution Ratio
N/A	N/A	N/A	N/A	N/A	N/A	N/A	1:100


  

Color Settings	Transparency (%)	Adjusted Hysteresis (%)
0	94	


# Transfer Data from iPad (Finder) – 1/2

Please refer to the “Transfer Data From iPad (iTunes)” if your computer is PC or MacOS Mojave or earlier.


1. Connect your iPad to your computer using the USB cable that comes with your device. Also, make sure that your MacOS Catalina is updated to the newest version.
2. Enter “Finder,” find your iPad in the left sidebar, and click.



3. Click the “Trust” button to connect your iPad mini.




4. You might see a prompt on the iOS device asking you to Trust This Computer. Press Trust to continue.




# Transfer Data from iPad (Finder) – 2/2

5. Click the “Files” label on the top.



6. Find the iSperm app on the application list, and select all the files under the application.



7. Drag the files into the target folder.




















# 06

## Frequently Asked Questions (FAQ)

# Can all types of frozen semen be analyzed?

Generally, there are three types of frozen semen.


Frozen sample with **unclear** extender **CANNOT** be analyzed by iSperm.

	Optically clear extender	Unclear extender	Sex-sorted semen
			
			
	USA 11HO11432	USA 1JE576	USA 511HO11314
iSperm			
Microscope			
CASA			
	(O)	(X)	(O)

# If Cover Chip is contaminated...

- If Cover Chip is contaminated, it is okay to use tissue paper to wipe the cover surface gently.
- It is recommended to keep the table surface, Base Chips, and Cover Chips clean.


## Microscopic view of dust/fiber on the Cover Chip.



# If bubbles are visible...

- **Re-sampling** is required if bubbles are visible. (See the picture below, Bubbles hinder the analysis significantly.)
- Bubble trapping is usually associated with the mixing process.
- Mix the semen gently can prevent it from forming bubbles.  
(It may take some practice to be skilled at bubble-free mixing.)

**Microscopic view of bubbles.**




# If image is blurry...

- Check if Base Chip and Collector are intact and locked correctly.
- Check if Sample Collector is locked correctly.
- Check if Cover Chip is clean.
- If the image is blurry in every chip, the lens may have been contaminated.

**Please get in touch with your distributor for technical service.**

## Microscopic view of blurry image.




# If sperm aren't distinguishable


## Possibilities:

- Sample could be too concentrated, and further dilution is required.
- Base Chip & Cover Chip aren't locked correctly to enclose a thin layer. It leads to multi-layers of sperm cells instead of a thin layer. In this case, re-sampling is required.

Single layer (O)



Multiple layers (X)




# Causes of Variation

## Possibilities:

- Mix Semen Uniformly


It is critical to make sure the sperm cells are distributed uniformly. See below examples of non-uniform semen. Further mixing is required.



**Mix semen thoroughly before Every sampling.**

- Time gap (delay) on Cover-Base Lock-in

It is important to **lock the Base Chip to the Cover Chip as soon as the sample is attached**. Experiments show that deviation exists with a 10-second time gap.



**After semen is attached, complete the sampling as soon as possible.**

# 07



## Validation








## Concentration Validation: iSperm vs. Nucleocounter & Hemocytometer





## Motility Validation: iSperm vs. CASA

The motility is compared at concentraton of 20-60M/ml.






## Concentration Validation: iSperm vs. Hemocytometer & CASA




## Motility Validation: iSperm vs. CASA

The motility is compared at concentration of 20-60M/ml.






## Concentration Validation: iSperm vs. CASA





## Motility Validation: iSperm vs. CASA

The motility is compared at concentraton of 20-60M/ml.






## Concentration Validation: iSperm vs. Hemocytometer & CASA




## Motility Validation: iSperm vs. CASA

The motility is compared at concentration of 20-60M/ml.






## Concentration Validation: iSperm vs. Hemocytometer & CASA




Pearson Correlation: 0.962  
Number of Samples: 110  
R<sup>2</sup>: 0.9843



Pearson Correlation: 0.967  
Number of Samples: 80  
R<sup>2</sup>: 0.9738

## Motility Validation: iSperm vs. CASA


The motility is compared at concentraton of 20-60M/ml.



Pearson Correlation: 0.973  
Number of Samples: 90  
R<sup>2</sup>: 0.9467




## Concentration Validation: iSperm vs. Hemocytometer



## Motility Validation: iSperm vs. CASA

The motility is compared at concentraton of 20-60M/ml.



# 08

## Certificates

## Certificate of Compliance



No. 0L160311.AB0558

Technical Construction File no. AID-2016001-A1

Certificate's Holder: Aidmics Biotechnology Co., Ltd.  
Rm. 1, 11F., No.171, Sec. 3, Roosevelt Rd., Da'an Dist., Taipei City 106, Taiwan (R.O.C.)

Certification ECM Mark:



Product: Sperm Analyzer  
Model(s): ADBISP

Verification to: Standard:  
EN 61000-6-2:2005/AC:2005,  
EN 61000-6-4:2007/A1:2011,  
EN 61000-4-2:2009, EN 61000-4-3:2010,  
EN 61000-4-8:2010

related to CE Directive(s):  
2014/30/EU (Electromagnetic Compatibility)

**Remark:** The product(s) has been verified on a voluntary basis. The product(s) satisfies the requirements of the Certification Mark of ECM, in reference to the above listed Standard(s). The above Compliance Mark can be affixed on the product(s) accordingly to the ECM regulation about its release and its use. The regulation can be found at [www.entecerma.it](http://www.entecerma.it). This Certificate of Compliance can be checked for validity at [www.entecerma.it](http://www.entecerma.it)

This verification doesn't imply assessment of the production of the product(s).

Additional information, clarification about the **CE** Marking:



We attest that a TCF for the **CE** Marking process is in place. Whereas the Manufacturer is Responsible to start the **CE Marking Certification Procedure** and to perform all the necessary activities, as required by the Directive before placing the **CE** Mark on the product(s).

Date of issue 11 March 2016

Chief Manager  
Tim Mahon



Expiry date 10 March 2021

Deputy Manager  
Viola Miller



Ente Certificazione Macchine Srl

Via Ca' Bella, 243 – Loc. Castello di Serravalle – 40053 Valsamoggia (BO) – ITALY  
☎ +39 051 6705141 ☎ +39 051 6705156 ✉ [info@entecerma.it](mailto:info@entecerma.it) 🌐 [www.entecerma.it](http://www.entecerma.it)



## CERTIFICATION

**Applicant** : Aidmics Biotechnology Co., Ltd.  
**Address** : Rm.1, 11F., No.171, Sec.3, Roosevelt Rd., Da'an Dist., Taipei City 106, Taiwan (R.O.C.)  
**Manufacturer** : Aidmics Biotechnology Co., Ltd.  
**Address** : Rm.1, 11F., No.171, Sec.3, Roosevelt Rd., Da'an Dist., Taipei City 106, Taiwan (R.O.C.)  
**Description of EUT** : Sperm Analyzer  
**Trade Name** : N/A  
**Model Number** : ADBISP  
**Product Series** : N/A  
**Type of Test** : FCC Part 15 Subpart B  
**Technical Standard** : **Emission**  
 FCC Part 15 : Subpart B Class A  
 CISPR 22 : 2008 Class A  
**Report Number** : HA160126-FD  
**Receipt Date** : 24-FEB-2016  
**Issue Date** : 08-MAR-2016  
**Test Result** : **Compliance**

The above equipment was tested by *HongAn TECHNOLOGY CO., LTD.*, for compliance with the requirement set forth in the FCC Rules and Regulation Part 15, Subpart B and the measurement procedures were based on ANSI C63.4.

**Note :**

1. The results of the test report relate only to the sample tested.
2. The test report shall not be reproduced without the written approval of *HongAn TECHNOLOGY CO., LTD.*

**Approved by:**

*Adam Yang*

Adam Yang / Section Manager



**HongAn TECHNOLOGY CO., LTD.**  
 NO.15-1, CWEISHUH KENG, CWEIPIN VILLAGE,  
 LINKOU DIST, NEW TAIPEI CITY, TAIWAN, R.O.C.

**TEL :** +886-2-26030362  
**FAX :** +886-2-26019259  
**E-mail :** haf@ms19.hinet.net

**BSMI Registration No.:** SL2-IN-E-0023, SL2-IS-E-0023,  
 SL2-A1-E-0023, SL2-R1-E-0023,  
 SL2-R2-E-0023, SL2-L1-E-0023

**FCC Designation No.:** TW1071, TW1163  
**TAF Accreditation No.:** 1163  
**VCCI Registration No.:** R-2156, C-2329, T-219, G-686







**Aidmics Biotechnology**

[www.aidmics.com](http://www.aidmics.com)  
[service@aidmics.com](mailto:service@aidmics.com)